

Flo-Assist® II Pressure and Flow Enhanced Irrigation



Control . . . Efficiency . . . Disposable

FLO-ASSIST® II provides control of irrigation flow during endoscopic procedures utilizing gravity fluid flow. A convenient hand or foot pedal gives the physician control over the timing and increase of irrigation pressures during the procedure.

Item Number: 7-510-49

**FLO-ASSIST® II Tubing Set
with disposable foot pedal**

Sold sterile, 10 sets per box at \$545.00

DISPOSABLE TUBING SET WITH REQUIRED CONNECTORS

- Allows physicians to control irrigation flow utilizing gravity fluid flow
- Completely disposable set helps control costs
- Simple, sterilized, disposable tubing set-up with ratchet clamp at bag end and roller clamp at spiked end
- Contains a drip chamber for visualization of fluid-flow
- Foot pedal allows for hands-free operation
- Compatible with a Foley catheter and CYSTO tubing set

FLUID MANAGEMENT CHAMBER

- Optional hand or foot-assisted control of irrigation flow
- Patented design includes a special fluid management chamber that automatically preloads with fluid
- Proprietary system keeps a surcharge of fluid while continuously delivering standard flow
- Bolus fluid delivery improves scope insertion and allows for a clear field of view

MODULAR DESIGN

- Configured for post-operative irrigation through a Foley catheter
- Compatible with various bag spike connectors for operating procedures using two or four spike systems
- Easy-to-implement and train product conversion in any surgical suite

To order the FLO-ASSIST® II Fluid Delivery System and Tubing Set for urologic and arthroscopic procedures contact NTI customer service at customerservice@ntisurgical.com or call 1-800-348-0424

Northgate Technologies Inc.
1591 Scottsdale Court
Elgin, IL 60123 USA

1-800-348-0424 • 01-847-608-8900 (Outside of the USA) • 01-847-608-9405 (FAX) • www.ntisurgical.com

Northgate Technologies, Inc. is a U.S.-based company™/® trademarks are registered trademarks of or are used under license by NORTHGATE Technologies, Inc.
© 2018 Northgate Technologies, Inc. All Rights Reserved.



Northgate Technologies, Inc.