



Convert your expense into manageable savings with the NEBULAE™ I System




NEBULAE™ I Insufflator

+



Insufflation Tubing Set

+



TAP Sensing Tubing Set

+



In-Line Gas Warmer

NEBULAE™ I

Pairs advanced technology with our intimate knowledge of surgeon and OR staff needs to provide the maximum value for everyday to advanced laparoscopic procedures.

- + Gynecology
- + Urology
- + Pediatric
- + General Surgery
- + Bariatric Surgery
- + Robotic Surgery

Commercial Offerings

RENTALS (M7-650-00)



RENT-TO-OWN (M7-650-00)


Option to pay difference on the machine after 6 months of rental

PURCHASE (7-650-00)

Rentals and Rent-to-own are only available on the NEBULAE™ I insufflator

Create a secure and robust pneumoperitoneum

- ✓ **Use standard trocars**
Patented high-flow luer connector is designed to work with any standard trocar with a luer connector.
- ✓ **Reusable CO₂ Warmer**
Delivers CO₂ gas at a near-body temperature.
- ✓ **5 different modes of operation**
General Laparoscopic, Bariatric, Pediatric, Vessel Harvest, Custom.




* Data on file. ** In use with TAP Sensing Tubing Set and in passive smoke evacuation setup.



Insufflation Tubing Set

- Provides 50 LPM of CO₂ gas flow all the way to the end of the tubing set
- Hydrophobic, 0.12µ ULPA filter with 99.999% efficiency
- Patented connector design ensures consistent performance and quality on both proximal and distal end during CO₂ delivery



TAP Sensing Tubing Set

- Continuously monitors pneumo pressure without interrupting CO₂ gas flow
- Stabilizes pneumo during any laparoscopic procedure that requires passive smoke evacuation at a low or high rate
- Minimizes the need to adjust the robotic arms and instruments due to the loss of pneumo



In-Line Warmer Cable

- Warms CO₂ insufflation gas near the point of entry
- Delivers CO₂ at near-body temperature
- Supports flow rates up to 50 LPM
- Use up to 100 times with proper cleaning and sterilization

NEBULAE™ I System

PRODUCT	DESCRIPTION	QTY.	PRODUCT NO.
NEBULAE™ I Laparoscopic Insufflator	50 L/minute, High Flow Insufflator	1	7-650-00 Purchase M7-650-00 Rental
Insufflation Tubing Set	Hydrophobic, 0.12µ ULPA filter with 99.999% efficiency Length: 10 ft.	10 / case Sterile	7-510-52
TAP Sensing Tubing Set	Length: 10 ft.	10 / case Sterile	7-510-31
In-Line Warmer Cable	In-Line Warmer, Heater Assembly Cable Length: 10 ft.	1	6-820-00

For more information on Northgate Technologies, Inc.:

Call **1-800-348-0424** inside U.S.A. | 01-847-608-8900 Fax: 01-847-608-9505 outside of the U.S.A
Visit [ntisurgical.com/products/insufflation](https://www.ntisurgical.com/products/insufflation) | E mail customerservice@ntisurgical.com

Northgate® Technologies, Inc. is a U.S.-based company™/® trademarks are registered trademarks of, or are used under license by NORTHGATE® Technologies, Inc.

© 2018 Northgate Technologies, Inc. All Rights Reserved.

82-00050-1 REV. B

