


NEW PRODUCT

Insufflation Systems
Accessory Product
from NORTECH®

CO₂ Shuttle Valve System

Avoid the stress and havoc of searching for a tank of CO₂ insufflation gas at a critical time during surgery.

The CO₂ shuttle valve is much more than a brass "T" fitting. It is a precision valve designed to sense the gas bottle with higher pressure and will switch over automatically when a new bottle of gas is opened. This means there is no need for a surgical procedure to be interrupted causing a potential dangerous situation or delay.



Ask for.....
NORTECH
Cat. # 7-510-50

- A. CO₂ gas bottle yokes
- B. 4 foot long pressure hose
- C. The Shuttle Valve connected to 2, 1 foot long pressure hoses.

NORTECH®

Since 1972.....

*The Pioneers of High Flow
Insufflation Technology*



Insufflation Product

Northgate Technologies Inc. • 600 Church Road • Elgin, Illinois 60123 847 608-8900 • Toll Free 800 348-0424

® / ™ Registered trademarks and trademarks of or used under license by Northgate Technologies Inc.